# iCUBE/MICROLAB/RAMPOWER

# OWNERS' MANUAL

Please take a moment to read this manual prior to operation!






182 Ridge Road, Suite D, Dayton, N.J 08810





## **CONTROLLERS**

# **iCUBE** 11111111111 © RAM PRODUCTS, INC.





## Set-Up

- Place in a convenient location avoiding any moisture or chemicals.
- •All components should be connected prior to turning on the main power switch (#3).
  •Make sure the power cord is securely attached in the rear (#1).
- \*Controllers utilize a state of the art microprocessor chip for programming of speed control, overload and safety and much more. This includes a memory feature to remember the last user settings before the main power is turned off.
- •Turn off the main power switch when not in use to avoid power and electric current to
- Always check the handpiece carefully for smooth operation prior to each use.
- •Do not use with any other brands of equipment as brushless units are not compatible with other brands.
- **1** Power Outlet: Main power cord receptacle. The power cord is plugged here. Warning: The control box is designed for one Input AC voltage only. Make sure that the available electrical power is compatible with the voltage input requirement of the control box. Refer to #2 as required.
- 2 110/230V: Factory preset switch is 110 volts but can accommodate 230 volts.
- Main Power Switch: Turns the control box on and off. Switch should be in the "off" position when not in use.
- Foot Switch Control Outlet: The foot pedal plugs in here, either the on/off or variable speed foot pedal.
- 5 Handpiece Connector: The handpiece cord plugs in here.
- 6 Speed Display: LED speed display and error alerts.
- 2 Speed Control Dial: 0-35,000 rpm or 0-45,000 rpm. Used to adjust or vary the desired speed. In hand mode, turn the dial clockwise to increase the speed and counterclockwise to reduce the speed. While using in foot pedal mode, the speed dial setting can also limit the top or maximum speed for greater control and safety.
- B Hand/Foot Button: When set in the hand position, the hand indicator light will turn on, the speed is completely controlled by the speed dial knob, when set in the foot pedal mode, the foot indicator light will turn on, the speed is controlled by the settings of the knob when using the on/off pedal and the speed is completely variable when using the variable foot pedal standby when flashing or blinking.
- Motor Button: Engages the handpiece for operation and use or can disengage for stopping or taking a rest from operation while remaining ready.
- Forward/Reverse Button: Controls handpiece direction to forward or reverse.

## **Specifications**

Model	Input	Output	MaxAmp	Rotate Speed	Power	L	W	Н	Wt.
iCUBE	AC110V/230V 50HZ-60HZ	DC0-32V	2A	0-35,000rpm/ 0-45,000rpm	50W	5"	4"	4 1/2"	3.8lbs
MICROLAB	AC110V/230V 50HZ-60HZ	DC0-32V	2A	0-35,000rpm/ 0-45,000rpm	60W	8"	4"	6"	4.4lbs
RAMPOWER	AC110V/230V 50HZ-60HZ	DC0-32V	2A	0-35,000rpm/ 0-45,000rpm	60W	6"	3 1/4"	5 1/4"	4.4lbs

### **Assembly**

- 1. Handpiece cord plugs into the lower middle socket(#5).
- 2. Foot Pedal plugs into the lower middle socket(#4).
- 3. Plug the power cord to an electrical outlet. Make sure that the outlet is properly grounded and compatible to the voltage input requirement:

NOTE: Please make sure that the power switch is in "off" position before plugging the power cord to an electrical outlet. The speed control must be set to minimum and check the voltage switch selection to be sure it is set at the proper voltage either 110 or 230 volts.

- 4. Turn the power switch (#3) to "on" position. The LED display (#6) should appear.
- 5. The handpiece is now ready for operation. Press Hand/Foot Button (#8)
  - a. HAND MODE Hand indicator light should glow. The speed of the motor is set by turning the control knob clockwise.
  - b. FOOT CONTROL MODE Foot indicator light should glow. •With the On/Off foot pedal, step on the foot pedal and gradually turn the speed control knob clockwise until the desired speed is attained. •With the variable foot pedal, gradually step on the pedal until the desired speed is attained.



## TECH2000/45000 HANDPIECE

- H37L1 Hand Cap Spacer 26, Roller Body Pin Hand Cap Body
- Chuck Housing C-Ring 5. Collar
- Bearing Fix Washer 6. 7. 1260 Wave Washer
- 8. 11260 Washer 9. Ball Bearing 1260ZZ
- 10. Bush
- 11. Ball Bearing 1480ZZ
- 12. Spring Sheet 13. Coil Spring
- 14. Test Bur
- 15. Collet Chuck 2.35
- 16. Chuck Housing 17. Chuck Spring
- 18. H37L Sealing O-Ring
- 19. Joint Shaft Pin 20. Joint Shaft
- 21. Connector Coil Spring
- 22. Connector
- 23. Bearing Bush Collar
- 24. Bearing Bush Collar Coil Spring

- H37L1 Hand Cap Cover 25, Roller Body
  - 27. Roller Guide 28. Set Ring
  - 29. Cam
  - 30. Slide Washer 31. Slide Body
  - 32. Ball 33. Motor Case C-Ring
  - 34. Motor Case Coil Spring 35. Motor Case Spring Sheet
  - 36 1030 Fix O-Ring 37. Motor Case Cover
  - 38. Motor Case
  - 39. Motor Case Bolt
  - 40. Magnet
  - 41. Copper Washer 42.1030 Wave Washer
  - 43. H37L Ball Bearing 1040ZZ 67. BR Brush Bolt 44,1030 Copper Washer
  - 45. 1030 Teflon Washer
  - 46.1030 Spacer 47. Shaft
  - 48. Core

- 49. Front Shaft Bush 50. Rear Shaft Bush
- 51. Commutator
- 52. Coil
- 53. Wedge
- 54. Fan 55. 830 Spacer
- 56. Teflon Washer
- 57. Bearing Housing 58. Back Bearing Sit
- 59. Brush Holder
- 60. Carbon Brush
- 61. Brush Spring 62. Brush Cap
- 63. Terminal(L)
- 64. Terminal(R) 65. Carbon Brush Bolt
- 66. Ball Bearing 1030ZZ
- 68. Motor Cord
- 69. Carbon Housing
- Cover (Black) 70. Spindle Ass'y
- 71. Armature Ass'v

## **OPTIMUS 35/OPTIMUS 45 HANDPIECE**

- Test Bur
- Collet chuck 2.35
- Nose Cap
- 4. Nose Cap O-Ring
- Chuck Housing C-Ring
- 6. Collar
- B/R Stop Collar 7.
- 814 Wave Washer
- 9. 814 Washer 2
- 10. 814 Washer 1
- 11. Ball Bearing (814ZZ)
- 12. Chuck Housing
- 13. Connection Pin
- 14. Spacer
- 15. Joint Shaft
- 16. Disc Spring
- 17. Spring Stopper
- 18. Spring Pin
- 19. O-Ring
- 20. Connector
- 21. Spring Spacer
- 22. Coil Spring1
- 23. Spring Sheet
- 24. Hand Cap Stopper
- 25. Hand Cap (Blue)

- Set Ring 27. Set Kev
- 28. Cam O-Ring
- 29. Cam
- 30. Ball
- 31. Moving Housing 32. Coil Spring 2
- 33. B/R Collar
- 34. Ball Bearing(1360ZZ)
- 35. B/R Stop Bolt 36. Motor Case C-Ring
- 37. Motor Case Coil Spring 38. Motor Case Spring Sheet
- 39. 1030 Fix O-Ring
- 40. Motor Case
- 41. Motor Case Cover (Blue)
- 42. Motor Case Bolt
- 43. Magnet
- 44. Magnet Guide 45. Copper Washer
- 46.1030 Wave Washer
- 47. Ball Bearing 1040ZZ 48, 1040 Copper Washer
- 49.1040 Teflon Washer 50, 1040 Spacer

- 51. Shaft Ø4 52. Core
- 53. Front Shaft Bush Ø4
- 54. Commutator
- 55. Coil
- 56. Wedge 57. Fan
- 58, 830 Spacer
- 59. Teflon Washer
- 60. Bearing Housing 61. Back Bearing Sit
- 62 Brush Holder
- 63. Carbon Brush
- 64. Carbon Brush Spring
- 65. Carbon Brush Cap 66. Terminal (L)
- 67. Terminal (R)
- 68. Carbon Brush Bolt
- 69. Ball Bearing 1030ZZ 70. Bearing Brush Bolt
- 71. Motor Cord 72. Carbon Housing
- Cover (Blue) 73. Spindle Ass'y
- 74. Armature Ass'y 75. BHS1 Ass'v


## **MAINTENANCE**

## **Safety Precautions!**

- · Do not force any burs into the collet chuck and always check after insertion if the bur is securely
- locked by pulling the bur. Always operate at a safe speed for the size and weight of the bur. Start at 0 rpm while slowly increasing the speed while securely gripping the handpiece. Operating at higher speeds often does not cut any better or faster and can increase wear on the bearings in the handpiece as well as the cutting tool.
- · Never use bent, rusted or damaged burs.
- Follow basic safe speed operation according to the chart below and also the safety speed
- ratings of the manufacturer of the burs. Never operate without a bur properly locked in the handpiece.
- After usage, remove any excess dust or debris.
- · Do not operate near water or moisture.
- · Avoid any chemicals and disinfectants.

## Safety Speed Table Guidelines

This chart for allowable rotation speed can be slightly different depending on bur materials or bur type.





Select the proper tools and speed for the job.



ball bearing.



Visually inspect before use. Ensure the handpiece operates smoothly and without noise or vibration. DO NOT USE IF THERE IS ANY UNUSUAL NOISES OR VIBRATION! If the handpiece is in need of service,



contact Ram Products. When installing any kind of bur to your handpiece, insert it properly to the final depth. Proper mounting or installation should avoid any extended distance from the chuck to reduce vibration but be sure the bur is not touching the nose tip of the handpiece with distance to spin freely.



5 Use proper speed that has not exceeded the maximum safety speed. Exceeding the maximum safety speed can be dangerous or can damage the handpiece.



6 When using the handpiece, wear safety goggles/glasses.



7 Never use the product for any purpose other than the original intended use.



8 We suggest to store in a cool and dry place when not being used and to avoid humidity and moisture if possible.



9 Never leave your handpiece with the power switch turned on, an accident may occur if it is over-heated.

3

4

## WARRANTY PERIOD

Warranty is valid within the United States only.

ELECTRICAL EQUIPMENT - These products are under warranty for a period of ONE YEAR from the date of original installation by the purchaser.

HANDPIECES - These products are under warranty for a period of SIX MONTHS from the date of original installation by the purchaser. IT DOES NOT COVER MISUSE, NORMAL BEARING, OR CARBON BRUSH WEAR.

IN HOUSE OR FACTORY REPAIR - Repairs made under warranty will be honored until the completion of the original product warranty period. SHIPPING AND HANDLING CHARGES ARE NOT COVERED WITHIN THE WARRANTY.

## PRODUCT WARRANTY

All returned items shall be free from defects/alterations in material, workmanship and original product design.

- The product shall have been used under the conditions for which it has been intended, and in a manner for which it has been designed.
- The product shall have been operated in a manner for which it was intended and in the proper environmental conditions.
- The customer must maintain the product as specified within the individual product instruction booklet.

Electrical and mechanical equipments and accessories which are mishandled, used incorrectly, and/or used for a purpose other than that for which they were designed and intended, will void the product warranty.

## WARRANTY CLAIMS

Should your product require service attention within the warranty period, we will elect to make the necessary repair or replacement at no charge and providing the product is delivered prepaid to the proper address. However it is recommended that returns be adequately insured against carrier loss or damage.

All claims should describe the specific problem for accurate analysis and evaluation, the return date of the warranty card along with the supply house and/or dealer name.

No product replacement is permitted without prior authorization.

PRODUCT:	
PURCHASE FROM:	
PURCHASE DATE:	
INSTALLATION DATE:	

## WARRANTY REGISTRATION

## REGISTER ONLINE at www.ramprodinc.com/warranty

In order to validate your warranty, this card must be filled out completely and mailed within 10 days of purchase/installation.

## PLEASE PRINT OR TYPE

CITY:	
PRODUCT:	ZIP CODE:
SERIAL NO(S).:	
DEALER OR SUPPLY HOUSE NAME	AND ADDRESS:
The short survey below is designed of product supplied to our customer complete the survey.	
Did the item arrive in a satisfactory of If NO please explain:	condition? YES or NO?
Did the instruction manual provide a installation and operation? YES or N	all details clearly for proper O? If NO please explain:
What factor(s) were involved in your O Advertisement O Distributor Recommendation	Previous Purchase of a Ram Product
Any feedback, comments or suggest	tions:

Thank you for your feedback.







