

Univox® CLS-1

Compact loop driver with Dual Action AGC

Features

- Compact & powerful induction loop amplifier for use with TV/digibox/ satellite/VHS/DVD/stereo/MP3/ microphone etc
- Univox Dual Action AGC for smooth quality sound
- Volume thumb wheel, bass and treble controls
- Rugged metal case with effective cooling for long lasting use
- Three inputs: Dedicated sensitive mic input (12V phantom power/signal level potentiometer) a mic/line input and a line input
- Easy connection to loop pad, sofa or room loop
- Separate headphones output
- Blue LED's for indication of power connection, signal source level and loop current

A compact powerful loop amplifier suitable for cashier offices, information counters and homes.

Univox® CLS-1 is a compact yet powerful loop amplifier for wireless listening through your hearing aid. CLS-1 is developed for maximum power, versatility and ease of use and is equipped with Univox automatic gain control (Dual Action AGC) for a smooth sound experience. The thumb wheel volume control, with digits for easily memorized setting of the right field strength level, along with bass and treble controls, secures pure sound quality. The metal case and effective cooling ensures long lasting problem free use.

Connect CLS-1 to your TV, stereo, MP3, microphone or any other sound source that you want to hear clearly. A super sensitive dedicated microphone input with 12V phantom power and input level potentiometer, a combined mic/line input and a line input makes it easy to connect your CLS-1 to virtually any sound source you like. You can monitor your CLS-1 with three indicators in the front panel with blue LED's for power connection, signal source level and loop current.

Choose between loop pad, sofa loop and room loop cable for use with your CLS-1 and cover exactly the area you want. Place the loop pad in your favourite chair, the sofa loop underneath your sofa or the room loop around the walls of your room and tune your hearing aid into T position. With the headphones output you can even listen to the magnetic field acoustically.

For check out counter and reception desk applications, please refer to CTC-120/121.









Sofa loop



Room loop with 50 clips



Headphones output

Area of coverage 50m² room loop according to IEC 60118-4

Loop output Max. current 9App, for loop pad (max. 2 loop pads), sofa loop cable or room loop cable

> Max. voltage 6.6Vpp (open output) 70-15,000Hz (±3dB) Frequency response

<1% Distortion

Connection 2 pin DIN (loop pad and sofa loop cable) or spring clip terminal (loop cable)

Short circuit protection 100% short circuit protection, constant current down to 0 0hm;

Internal excessive current and temperature protection with auto-reset

Switches automatically to voltage feedback for correct frequency response,

70-15,000Hz (±3dB), when using headphones

Connection 3.5mm jack

Mic input 2-170mV/5k Ω (phantom voltage 12V). Separate high gain amplifier with AGC Inputs

> Connection 3.5mm jack

Mic/Line & Line inputs 10mV-4.3V/5kΩ (phantom voltage 12V in Mic/Line input) knee-point setting

Connections 3.5mm jack

Dual action AGC Dynamic range >70dB

> Attack / Release time 2-500ms / 0.5-20dB/s

Trimmers Bass control 0 - +12dB, potentiometer at front panel

Treble control 0 - +9dB, potentiometer at front panel Volume control Graded 0-9, thumb wheel at front panel

Optical indication Mains connection Blue LED

Loop current Blue LED Blue LED Input signal level

Other information **Dimensions** 150x92x25mm (LxWxH)

Weight 395g including power supply

Part No 202052 (EU) / 202054 (UK) / 202056 (US) / 202058 (AUS)

Note

Power supply Technical data may vary due to different power supplies

Dual action AGC Variable attack/release times holding peak and average level constant

Accessories

Part No Description 900734 Switched power supply, 230/15V 18VA (EU) 900738 Switched power supply, 230/15V 18VA (UK) 900740 Switched power supply, 120/15V 18VA (US) 900744 Switched power supply, 230/15V 18VA (AUS) 283002 Loop pad, 400x400x5mm, 2 pin DIN, grey, cable=10m 283210 Extension cable for loop pad, 2 pin DIN, black, cable=10m

283220 Sofa loop cable incl. ext. lead, 2 turns, cable connector, 2 pin DIN, black, cable=8+10m

283240 Room loop cable roll, 0.75 mm², incl cable clips, light grey, cable=30m

Adaptor cable for car cigarette lighter for Univox CLS-1 289002 281015 Audio cable, 3.5mm male to two RCA male, cable=1.5m

281041 HDA-2M, D/A converter incl. optical cable

282007 Ramification cable, 2 pin DIN male to two 2 pin DIN female, cable=0.2m

13A, Electret mic, velcro adhesive, 3.5mm male, cable=3m 242401 241501 13T, Electret microphone, condenser capsule, boundary

Holder for wall mounting, black aluminium

Kits

289015

Choice of Loop pad Sofa loop Room loop 202052P 202052S 202052R EU power supply UK power supply 202054P 2020545 202054R US power supply 202056P 202056S 202056R AUS power supply 202058P 202058S 202058R

202060C Car kit CLS-1, adaptor cable, 13A velcro adhesive microphone, tie-clip microphone and loop pad





