

Printed: Mon Nov 22 11:59:49 EST 2010



Viewing components of :21700

Item 7 on the 21700 is black

Item 7 on the 21760 is green

Part Number	Ref. #	Qty	Description
<b>03100-LED</b>		1	Led Lamp Upgrade Kit 3.5V Diag
<b>03100-U</b>	1	6	3.5V Halogen Lamp
<b>106008</b>	2	1	Screw 3-56 .100 Prism Flat Slotted
<b>203009</b>	3	1	Washer .240/.154 Spring
<b>203010</b>	4	2	Lens Holder Washer
<b>209012</b>	5	1	Screw 3-56 .182 Pan Slotted
<b>209025-1</b>	6	1	Spec Holder
<b>209026-501</b>	7	1	Lens Holder Assy-black
<b>209026-502</b>	7	1	Lens Holder Assy-green
<b>217009</b>	8	1	Nut 15/32-32 .577 Locking Thumbscrew
<b>217010-1</b>	9	2	Lock Ring
<b>217023</b>	10	1	Body Sub-assy W/lite Pipe
<b>250004</b>	11	1	Lamp Collar

