



# **Ultra Comfort Air Spring Stool**

### Model 6749

All UMF Medical stools are available in 16 best-in-class seamless upholstery colors with Prefixx® Protective Finish to assist with infection prevention and for full exam room cohesion.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



### **Colors**



## Key features at a glance

- Ultra Comfort seamless seat with anti-microbial guard
- > 5 leg black plastic base
- 2" (5.1cm) twin wheel casters
- Air spring height adjustment
- Available in all UMF upholstery colors
- 3-year limited warranty

### **Specifications**

Seat diameter: 16" (40.6cm) Min. seat height: 17.25" (43.8cm) Max. seat height: 22.5" (57.2cm) Net Weight: 13lb (5.9kg) Shipping Weight: 18lb (8.2kg) Shipping Dimensions: L23xW24xH8 (58.4cm x 61cm x 20.3cm)

#### Accessories

6749F – Fire retardant foam seat cushion 150 – Stool glides — set of 5 111 – Auto-lock occupied casters

